



Single-use Containers for



HUMAN VENOUS BLOOD SPECIMEN COLLECTION

SHANGHAI KINDLY(KDL) GROUP
NANCHANG KINDLY MEDICAL TECHNOLOGY CO.,LTD



Specimen type and preparation method:

Lid Color	Clinic Application	Specimen Type	Steps of Specimen Preparation	Additive
	Biochemical test, blood bank test	Blood serum	After blood collection, you don't need to invert and mix, let it stand for 1 hour and centrifuge	Nothing
	Serum test	Blood serum	Recommended minimum time is 15-30 minutes 1500~2200g(RCF)for 10 minutes at 18-25°C	Coagulant
	PCR check	Blood plasma	Immediately after blood collection invert and mix 8 -10 times, centrifuge	Separation gel /EDTA.K2
	Blood conventional check, full blood test	Whole blood	Immediately after blood collection, invert and mix 8-10 times, mix sample before test	EDTA.K2 or EDTA.K3
	Quick biochemical test of blood, serum biochemical test	Blood serum	Recommended coagulation time is 15-30 minutes 1500~2200g(RCF)for 10 minutes at 18-25°C	Clot Activator & Gel
	Quick biochemical test of blood plasma, hemorrheology test	Blood plasma	Immediately after blood collection. invert and mix 8-10 times, centrifuge	Heparin sodium or heparin lithium
	Test of blood solidification, PT/APTT, test of solidification blood gene	Blood plasma	Immediately after blood collection. invert and mix 8 times,centrifuge	Na-Citrate:3.2% blood sample=1:9
	Blood sugar test	Blood plasma	Immediately after blood collection. invert and mix 8-10 times,centrifuge ≤1300g(RCF) for 10 minutes at 18-25°C	Sodium Fluoride Sodium Heparin
	Erythrocyte sediment rate test	Whole blood	Gently inverted 180in 8-10 times	Na-Citrate:3.8% blood sample=1:4
	Blood bank research, HLA typing, DNA paternity testing and other research	Whole blood	Gently mixed 180 degrees upside down for 8-10 times ≤1300g(RCF) for 10 minutes at 25°C	Anti-coagulant ACD(Acid Citrate Dextrose)
	Separation of undiluted plasma,molecular diagnostic test	Blood plasma	Immediately after blood collection, gently reverse blood vessel mixing 5-10 times	K2 EDTA & Gel
	Quick biochemical test of blood plasma,biochemical test	Blood plasma	Immediately after blood collection. invert and mix 8-10 times,centrifuge 1500-2000g(RCF) for 10 minutes at 18-25°C	Lithium Heparin & Gel

KDL Serum Tubes																																	
KDL serum tubes could produce(no additive) high quality serum, it doesn't contain clot activator with a gel barrier.																																	
Specification: No Additive.																																	
Cap Colour: Red & White.																																	
<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr style="background-color: #f2f2f2;"> <th>Ref.</th> <th>Volume</th> <th>Size(mm)</th> <th>Material</th> </tr> </thead> <tbody> <tr><td>6656</td><td>3ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6657</td><td>5ml</td><td>13x100</td><td>Glass</td></tr> <tr><td>6658</td><td>6ml</td><td>13x100</td><td>Glass</td></tr> <tr><td>6001</td><td>2ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6002</td><td>3ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6005</td><td>5ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6006</td><td>6ml</td><td>13x100</td><td>PET</td></tr> </tbody> </table>	Ref.	Volume	Size(mm)	Material	6656	3ml	13x75	Glass	6657	5ml	13x100	Glass	6658	6ml	13x100	Glass	6001	2ml	13x75	PET	6002	3ml	13x75	PET	6005	5ml	13x100	PET	6006	6ml	13x100	PET	<p>Clotting times:</p> <p>The recommended minimum coagulation time for serum tubes from patients who have not been treated with anticoagulant therapy is 60 minutes</p>
Ref.	Volume	Size(mm)	Material																														
6656	3ml	13x75	Glass																														
6657	5ml	13x100	Glass																														
6658	6ml	13x100	Glass																														
6001	2ml	13x75	PET																														
6002	3ml	13x75	PET																														
6005	5ml	13x100	PET																														
6006	6ml	13x100	PET																														



KDL Serum Tubes																													
To obtain serum samples from plastic tubes, coagulation activator needs to be added in the tube, since the plastic surface alone is insufficient to activate the blood coagulation within certain time, KDL plastic serum tubes are processed with silicon on the inner surface, which will accelerate the coagulation.																													
Specification: Serum with clot activator(silicaparticles)																													
Cap Colour: Red & Orange																													
<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr style="background-color: #f2f2f2;"> <th>Ref.</th> <th>Volume</th> <th>Size(mm)</th> <th>Material</th> </tr> </thead> <tbody> <tr><td>6014</td><td>3ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6018</td><td>5ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6019</td><td>6ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6662</td><td>3ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6664</td><td>5ml</td><td>13x100</td><td>Glass</td></tr> <tr><td>6665</td><td>6ml</td><td>13x100</td><td>Glass</td></tr> </tbody> </table>	Ref.	Volume	Size(mm)	Material	6014	3ml	13x75	PET	6018	5ml	13x100	PET	6019	6ml	13x100	PET	6662	3ml	13x75	Glass	6664	5ml	13x100	Glass	6665	6ml	13x100	Glass	<p>Clotting times:</p> <p>The recommended minimum time for the coagulation of serum tubes from patients who have not been treated with anticoagulants is 15-30 minutes</p> <p>Centrifuging conditions:</p> <p>1500~2200g(RCF) for 10 minutes at 18-25°</p>
Ref.	Volume	Size(mm)	Material																										
6014	3ml	13x75	PET																										
6018	5ml	13x100	PET																										
6019	6ml	13x100	PET																										
6662	3ml	13x75	Glass																										
6664	5ml	13x100	Glass																										
6665	6ml	13x100	Glass																										



KDL K2 & K3 EDTA Tubes																																																																									
EDTA salts (ethylenediamine tetracetic acid) are used in hematological studies to anticoagulate whole blood because the cellular components of blood could be well preserved by EDTA particularly. It is an anticoagulant because it forms a complex with metal ions such as calcium that inhibits the clotting cascade.																																																																									
EDTA anticoagulation is irreversible. When filling levels are correct, the concentration of EDTA in KDL tubes is 1.8mg/ml of whole blood, as recommended by the ICSH (International Council of Hematological Societies). The ICSH recommends the use of EDTA dipotassium (K2EDTA) for hematological examination. KDL plastic tubes are available with spray drying K2EDTA and K3EDTA*																																																																									
Specification: K2 EDTA	Specification: K3 EDTA																																																																								
Cap Colour: Purple	Cap Colour: Purple																																																																								
<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr style="background-color: #f2f2f2;"> <th>Ref.</th> <th>Volume</th> <th>Size(mm)</th> <th>Material</th> </tr> </thead> <tbody> <tr><td>6049</td><td>2ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6051</td><td>3ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6053</td><td>4ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6055</td><td>5ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6057</td><td>6ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6676</td><td>2ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6677</td><td>3ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6678</td><td>4ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6679</td><td>5ml</td><td>13x100</td><td>Glass</td></tr> <tr><td>6680</td><td>6ml</td><td>13x100</td><td>Glass</td></tr> </tbody> </table>	Ref.	Volume	Size(mm)	Material	6049	2ml	13x75	PET	6051	3ml	13x75	PET	6053	4ml	13x75	PET	6055	5ml	13x100	PET	6057	6ml	13x100	PET	6676	2ml	13x75	Glass	6677	3ml	13x75	Glass	6678	4ml	13x75	Glass	6679	5ml	13x100	Glass	6680	6ml	13x100	Glass	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr style="background-color: #f2f2f2;"> <th>Ref.</th> <th>Volume</th> <th>Size(mm)</th> <th>Material</th> </tr> </thead> <tbody> <tr><td>6066</td><td>2ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6068</td><td>3ml</td><td>13x75</td><td>PET</td></tr> <tr><td>6072</td><td>5ml</td><td>13x100</td><td>PET</td></tr> <tr><td>6683</td><td>2ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6684</td><td>3ml</td><td>13x75</td><td>Glass</td></tr> <tr><td>6685</td><td>5ml</td><td>13x100</td><td>Glass</td></tr> </tbody> </table>	Ref.	Volume	Size(mm)	Material	6066	2ml	13x75	PET	6068	3ml	13x75	PET	6072	5ml	13x100	PET	6683	2ml	13x75	Glass	6684	3ml	13x75	Glass	6685	5ml	13x100	Glass
Ref.	Volume	Size(mm)	Material																																																																						
6049	2ml	13x75	PET																																																																						
6051	3ml	13x75	PET																																																																						
6053	4ml	13x75	PET																																																																						
6055	5ml	13x100	PET																																																																						
6057	6ml	13x100	PET																																																																						
6676	2ml	13x75	Glass																																																																						
6677	3ml	13x75	Glass																																																																						
6678	4ml	13x75	Glass																																																																						
6679	5ml	13x100	Glass																																																																						
6680	6ml	13x100	Glass																																																																						
Ref.	Volume	Size(mm)	Material																																																																						
6066	2ml	13x75	PET																																																																						
6068	3ml	13x75	PET																																																																						
6072	5ml	13x100	PET																																																																						
6683	2ml	13x75	Glass																																																																						
6684	3ml	13x75	Glass																																																																						
6685	5ml	13x100	Glass																																																																						
<p>*International Council for Standardization in Haematology and Expert Panel on Cytometry, Recommendations of the International Council for the Standardization in Haematology for Ethylenediaminetetraacetic Acid Anticoagulation of Blood for Blood Cell Counting and Sizing. Am J Clin Pathol 1993;100: 371-372.</p>																																																																									



Tube mixing:
Correct mixing(8-10 inversions)of EDTA tubes immediately after blood collection is very important to avoid blood micro coagulation effectively.

KDL Advance tubes with Gel & Clot Activator

An inert gel is added in KDL advanced gel & clot activator tubes. The gel could separate the serum and the blood clot, which prevents the contamination of the serum from the separated cellular components during centrifugation.

For example, the serum for certain analytes such as potassium, phosphorus and glucose must be separated from the cells within a few hours - otherwise the results will be distorted considerably. Using KDL Advance tubes in routine clinical chemistry, the analytes such as potassium and glucose are still stable after a week of storage at 2-8°C.

Specification: Clot Activator & Gel

Cap Colour: Yellow

Ref.	Volume	Size(mm)	Material
6026	3ml	13x75	PET
6028	4ml	13x100	PET
6029	5ml	13x100	PET
6668	3ml	13x75	Glass
6669	4ml	13x100	Glass
6670	5ml	13x100	Glass

The main advantages of gel tubes:

- The barrier between serum and clotting is stable, therefore the analyte is more stable, the sample quality is better.
- Optimized workflow; Short centrifugation time, sample processing and archiving in the primary tube.
- Eliminate the possibility of misidentification due to the use of secondary tubes.

Clotting times:

The recommended coagulation time is 15-30 minutes.



Centrifuging conditions:

1500~2200g(RCF) for 10 minutes at 18-25°C

KDL Plasma Tubes

KDL plasma tubes are available with spray-dried sodium heparin or lithium heparin additives, these tubes can be used for clinical chemistry analysis. The heparin works as anticoagulant by developing an antithrombin complex, which inhibits thrombin and the activated factor X and thus prevents coagulation.

KDL plasma tubes with lithium heparin is spray dried on the surface of tube inner walls by special procedure, so that the additive is evenly distributed to achieve the best possible solubility. For clinical chemistry, lithium heparin is generally preferred over sodium heparin.

Specification: Sodium Heparin

Cap Colour: Green

Ref.	Volume	Size(mm)	Material
6102	3ml	13x75	PET
6104	4ml	13x75	PET
6106	5ml	13x100	PET
6108	6ml	13x100	PET
6696	3ml	13x75	Glass
6697	4ml	13x75	Glass
6698	5ml	13x100	Glass
6699	6ml	13x100	Glass

Specification: Lithium Heparin

Cap Colour: Green

Ref.	Volume	Size(mm)	Material
6122	3ml	13x75	PET
6124	4ml	13x75	PET
6126	5ml	13x100	PET
6128	6ml	13x100	PET
6702	3ml	13x75	Glass
6703	4ml	13x75	Glass
6704	5ml	13x100	Glass
6705	6ml	13x100	Glass



Tube mixing:

KDL heparin tubes should be mixed correctly (8-10 times inversions) and immediately after the blood sample has been taken, this is extremely important to avoid micro-clotting

KDL Plasma Tubes

KDL plasma tubes with separating gel for clinical chemistry are available with spray-dried lithium heparin additives. An inert gel could separate the serum and the blood clot preventing the contamination of the plasma from the separated cellular components during the centrifugation of the KDL Lithium Heparin & Gel tubes. The plasma for certain analytes such as potassium, phosphorus and glucose must be separated from the cells within a few hours, or else, the results will be significantly distorted. Using KDL Lithium Heparin & Gel tubes in clinical chemistry routine analytes, such as potassium and glucose are still stable after a week of storage at 2-8°C

Specification: Lithium Heparin & Gel

Cap Colour: Light Green

Ref.	Volume	Size(mm)	Material
6177	3ml	13x75	PET
6179	4ml	13x100	PET
6181	5ml	13x100	PET

The main advantages of gel tubes:

- The barrier between plasma and clotting is stable, therefore the analyte is more stable, the sample quality is better.
- Optimized workflow; Short centrifugation time, sample processing and archiving in the primary tube.
- Eliminate the possibility of misidentification due to the use of secondary tubes.

Effects of temperature:

KDL Lithium Heparin Gel tubes should be stored at 4-25°C and protected from direct sunlight during storage. Under centrifugation or low temperature conditions, the fluidity of the blood sample in the tube is reduced, The optimum temperature for separation of sediment and plasma is 20-25°C.



Tube mixing:

KDL Plasma tubes with separating gel should be mixed correctly (8-10 times inversions) and immediately after the blood sample has been taken, which avoids micro-clotting

Centrifuging conditions:

1500-2000g(RCF) for 10 minutes at 18-25°C

KDL Coagulation Tubes

KDL coagulation tubes contain a buffered sodium citrate solution and are used for examination of coagulation parameters, The tubes are available with a citrate concentration of 3.2%

The mixing ratio is 1 part citrate solution to 9 parts blood. Correct amount of blood collection is critical for correct coagulation analysis.

Specification: 3.2% Sodium Citrate

Cap Colour: Blue

Ref.	Volume	Size(mm)	Material
6082	2ml	13x75	PET/PP
6083	3ml	13x75	PET/PP
6688	2ml	13x75	Glass
6689	3ml	13x75	Glass

Coagulation conditions:

- Platelet-rich plasma: 1500-2000g(RCF) for 5 minutes at 18-25°C
- Platelet-poor plasma: Plastic tubes: 2000-2500g(RCF) for 10-15 minutes at 18-25°C
- Platelet-free plasma:> 3000g(RCF) for 15-30 minutes at 18-25°C

KDL recommends that glass tubes are not centrifuged at more than 2200g(RCF) in a swing-out rotor (for fixed angle rotor not more than 1300g(RCF))

Tube mixing:

Immediately after blood collection.invert and mix 8 times, centrifuge



KDL Plasma Tubes

KDL blood glucose analysis tubes can be used for clinical glucose analysis by adding sodium fluoride, heparin sodium, K2 EDTA and potassium oxalate. Glucose values of unpreserved blood samples drop rapidly after collection because glucose is metabolized by blood cells. The additive(sodium fluoride or sodium heparin) contained in the KDL blood glucosetube can stop enzymatic activity in the glycolytic pathway.

Specification: Sodium Fluoride
Sodium Heparin

Cap Colour: Gray

Ref.	Volume	Size(mm)	Material
6038	2ml	13x75	PET
6040	3ml	13x75	PET
6042	4ml	13x100	PET
6044	5ml	13x100	PET
6672	2ml	13x75	Glass
6673	3ml	13x75	Glass
6674	4ml	13x100	Glass
6675	5ml	13x100	Glass

Specification: Sodium Fluoride
Potassium Oxalate

Cap Colour: Gray

Ref.	Volume	Size(mm)	Material
6270	3ml	13x75	PET
6272	4ml	13x75	PET
6274	5ml	13x100	PET
6276	6ml	13x100	PET

Specification: Sodium Fluoride/K2 EDTA

Cap Colour: Gray

Ref.	Volume	Size(mm)	Material
6281	2ml	13x75	PET
6283	3ml	13x75	PET
6285	4ml	13x100	PET
6287	5ml	13x100	PET



Tube mixing:

KDL glucose tubes should be mixed correctly (8-10times inversions) and immediately after the blood sample has been taken

Centrifuging conditions:

≤1300g(RCF) for 10 minutes at 18-25°C

KDL ESR Tubes

KDL ESR tubes contain a buffered sodium citrate solution with concentration of 3.8%. The mixing ratio is 1 part citrate solution to 4 parts blood. Correct amount of blood collection is critical for correct Erythrocyte Sedimentation Rate analysis.

Specification: 3.8% Sodium Citrate

Cap Colour: Black

Ref.	Volume	Size(mm)	Material
6093	2ml	13x75	PET/PP
6094	3ml	13x75	PET/PP
6693	2ml	13x75	Glass
6694	3ml	13x75	Glass
6692	1.6ml	8x120	Glass
6695	5ml	13x100	Glass

Tube mixing:

ESR(Erythrocyte Sedimentation Rate)tubes should be gently inverted 180°in 8-10 times.



KDL ACD Tubes

KDL ACD tubes are suitable for blood bank research,HLA typing, DNA paternity testing and other research.It can effectively reduce the pre-analysis error of the specimen, improve the test accuracy,provide safety protection forpatients and medical staff,reduce nosocomial infection,The anti-coagulant ACD (Acid Citrate Dextrose) is used for the conservation of erythrocytes.

ACD is available in two solutions, A and B, each with different mixture ratios.

Studies have shown that the DNA preservation effect of ACD solution B is better than that of EDTA and heparin at room temperature.Under the protection of ACD solution, three cycles of freezing and thawing of samples had little effect on DNA.

Anti coagulant	ACD solution A	ACD solution B
Citrate	3.30mg/ml	1.89mg/ml
Citrate acid	1.20mg/ml	0.69mg/ml
Dextrose	3.68mg/ml	2.1mg/ml
Potassium Sorbate	0.03mg/ml	0.03mg/ml

The figures represent the final concentration in the blood in each case.

Specification: ACD Solution A

Specification: ACD Solution B

Cap Colour: Yellow

Cap Colour: Yellow

Ref.	Volume	Size(mm)	Material
6159	4ml	13x100	PET
6161	6ml	16x100	PET
6162	8ml	16x100	PET

Ref.	Volume	Size(mm)	Material
6167	4ml	13x100	PET
6169	6ml	16x100	PET
6170	8ml	16x100	PET



Tube mixing:

Gently mixed 180 degrees upside down for 8-10 times

Centrifuging conditions:

≤1300g(RCF) for 10 minutes at 25°C

KDL Plasma Preparation Tubes

KDL Plasma Preparation tubes are used for the separation of undiluted plasma from whole blood for molecular diagnostic test methods.

KDL Plasma Preparation tubes can ensure:

- Infectious samples could be handled safely
- The user will not be exposed to a blood sample in a test tube
- Plasma is prepared in closed KDL tubes and can be transported directly without the need for sub-loading and re-labeling from primary to secondary containers.
- Plasma quality is maintained, The gel barrier will prevent plasma from contacting with red blood cells to maintain plasma stability.

Specification: K2 EDTA& Gel

The viral load will remain stable for:

Cap Colour: White

- 6 hours - whole blood at room temperature
- 24 hours - separated plasma at room temperature
- 5 days - separated plasmarefrigerated at 4 °C

Ref.	Volume	Size(mm)	Material
6142	3ml	13x75	PET
6144	4ml	13x100	PET
6146	5ml	13x100	PET

Tube mixing:

Immediately after blood collection,gently reverse blood vessel mixing 5-10 times

Centrifuging conditions:

1500g(RCF)for 10 minutes at 18-25°C



SHANGHAI KINDLY(KDL) GROUP
NANCHANG KINDLY MEDICAL TECHNOLOGY CO.,LTD

F/7, Tower A, No.77 North Aixihu Road,Nanchang
High-Tech Industry Development Zone,Nanchang,
Jiangxi 330096.China

F: +86-791 86861216
T: +86-791 86861216
W: www.kdInc.com