



Integrated Aesthetic & Personal Care Solutions
OEM&ODM Manufacturer

Professional team | Patent technology | Medical standards



Shanghai Kindly Group

Nanchang Kindly Medical Technology Co.,Ltd.

Add: F/7,Tower A,No.77 North Aixihu Road, Nanchang High-Tec
Industry Development Zone, Nanchang, Jiangxi. China
Post Code:330096 Tel: +86-0791-86861216
Mob: +86-13922143599 Web: www.kdInc.com
E-mail: inquiry@kdInc.com

Shenzhen ENMIND Technology Co.,Ltd.

Add: 5th Floor, Block A, Defengsheng Building, No.41 Dabao
Road, Bao'an District 23, Shenzhen 518101, P.R.China.
Tel: 0755-2331 6007 Fax: 0755-2227 6344

NIR Skin Rejuvenation Beauty Device

Skin Brightening
and Rejuvenation

Reduce wrinkles

Control oil and
suppress acne

Lift and tighten
the skin



Technical Advantages

NIR Near-Infrared Light Technology
Multi-Frequency Microcurrent Technology
Photoelectric Integration Technology
6cm² Light Emission Area for rapid full-face coverage,
enhancing skincare efficiency.
Lamp tube lifespan of 3000 hours, usable for
up to 24 years with daily use once!

Product Information

Model Number: NIR-001

Wavelength Range: 900nm-1800nm

Product Input Power: 36W

Product Levels: 3 Levels

Mode: NIR Mode/EMS Mode/DER Mode

NIR + EMS+DER
Dual Efficacy with
Dual Electrodes for
More Effective Results

Technical Advantages

NIR Near-Infrared Light Technology
Multi-Frequency Multi-Electrode Microcurrent Technology
Photoelectric Integration Technology
6cm² Light Emission Area for rapid full-face coverage,
enhancing skincare efficiency. Built-in battery for greater portability
Lamp tube lifespan of 3000 hours, usable for
up to 24 years with daily use once!

Product Information

Model Number: NIR-002

Wavelength Range: 900nm-1800nm

Product Input Power: 36W

Product Levels: 3 Levels

Mode: NIR Mode/EMS Mode/DER Mode

Optimized Wavelength
900-1800nm
More Gentle and Effective
Light Therapy Experience



Supports ODM/OEM customized services.