



## Blood Collection Needle With Safety Device (V Type)

### Brief description of the product

The Blood Collection Needle with Safety Device are assembled by hub, needle tube, protect cap, safety device and rubber cap. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free. The Blood Collecting Needles intended to be used to blood or plasm collection.

### Characteristic


- ◆ Latex free
- ◆ The design of safe protector cap can shield the cannula after use, prevent needlepoint from stabbing nurse, patient or other people, which can effectively reduce infection.

### Detailed specifications

SIZE		COLOR	LENGTH			ID		Bevel angle		REF
OD (mm)	GAUGE	COLOR & CODE	3/4	1	1 1/2	RW	TW	Long bevel	Short bevel	
0.60	23G	blue								01-413
0.70	22G	black								01-414
0.80	21G	green								01-415

Remark: normal size in a above table, we can also make special size as clients' request.

### Instruction for use



### Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

### Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

### Main material

ABS,rubber, SUS304 Stainless Steel Cannula,Silicone oil

