



Cosmetic Needles

Brief description of the product

The Cosmetic Needles are assembled by hub, needle tube, protect cap. All materials meet medical requirement; sterilized by ETO, pyrogen-free.

The Cosmetic Needles is used for special injection tasks such as injecting filling material in Plastic surgery.

Characteristic

- ◆ Made of high quality stainless steel
- ◆ The needle tip smooth
- ◆ Size distinguished by the needle-hub's color for the clear recognition
- ◆ Production of all kinds of special needles required by our customers
- ◆ Sterile

Detailed specifications

Specifications	REF
1.2×70 TW	01-501
0.7×70 TW	01-502
0.7×50 TW	01-503
0.7×40 TW	01-504
0.6×50 TW	01-505
0.5×50 TW	01-506
0.5×40 TW	01-507
0.4×40 TW	01-508
0.3×25 RW	01-509

Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

ABS,rubber, SUS304 Stainless Steel Cannula,Silicone oil

