



Dental Needles

Brief description of the product

The Dental Needles are assembled by hub, needle tube, protect cap. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free. It's intended to be used anaesthesia in dentistry

Characteristic

- ◆ High quality stainless steel
- ◆ Custom needle style per specification

Detailed specifications

SIZE	COLOR	LENGTH		ID	BEVEL ANGLE		WORMTYPE		REF		
		Short	Long		Long bevel	Short bevel	METRIC	INCH			
OD (mm)	GAUGE	COLOR & CODE	1/2	5/8	1	1 1/4	1 1/2	RW	TW		
0.3	30G	purple									
		white	white							white	
0.40	27G	blue									
		yellow									
0.50	25G	yellow									
		orange									

Remark: normal size in above table, we can also make special size as clients' request.

Product Conformance

Conform to ISO 7885

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula, Silicone oil



1 Metric thread

2 Inch screw threads