



Vet Needles

Brief description of the product

The Vet Needles are assembled by hub, needle tube, protective cap. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free.

The Vet Needles are used together with Syringes or Infusion/Transfusion sets for medicine infusion and blood collection.

Characteristic

- ◆ Designed for veterinary use
- ◆ Sterile, non toxic, pyrogen free
- ◆ Autoclavable
- ◆ Needles and boxes are color-coded for easy identification
- ◆ Each needle is individually packaged and sealed
- ◆ Needle hub: copper hub, aluminum hub

Detailed specifications

SIZE	LENGTH	ID	BEVEL ANGLE		REF									
			RW	TW										
OD (mm)	GAUGE	1/3"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"	RW	TW	Long bevel	short bevel	
0.40	27G													01-1101
0.45	26G													01-1102
0.50	25G													01-1103
0.55	24G													01-1104
0.60	23G													01-1105
0.70	22G													01-1106
0.80	21G													01-1107
0.90	20G													01-1108
1.10	19G													01-1109
1.20	18G													01-1110
1.60	16G												white	01-1111
1.80	15G													01-1112
2.10	14G													01-1113

Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

SUS304 Stainless Steel Cannula, Silicone oil



Copper Hub



Aluminum Hub

