



Insulin Pen Needles

Brief description of the product

Insulin Pen Needles are assembled by needle hub, needle, small protective cap, big protective cap and other integral parts, for use with pre-diabetic insulin liquid filled insulin pen (such as the Novo Pen) for insulin injection. All these parts and materials used in this product can meet medical requirement, after ETO sterilization, pyrogen-free.

Characteristic

- ◆ Optimum treatment with maximum simplicity
- ◆ Sterile, non-toxic, non-pyrogenic
- ◆ Double point system
- ◆ Patient point with triple beveling for easy and atraumatic
- ◆ Cartridge point sharpened for a perfect punch with no fragmentation
- ◆ Compatible with all available pen models
- ◆ Easy flow of medication without any impurities

Detailed specifications

| Description | | | | | Catalog NO. |
|-------------|------|-----|------|--|-------------|
| 29G | 0.33 | 1/2 | 12mm | | 01-1201 |
| 30G | 0.3 | 1/4 | 4mm | | 01-1202 |
| 31G | 0.25 | 1/4 | 6mm | | 01-1203 |
| 31G | 0.25 | 1/3 | 8mm | | 01-1204 |
| 32G | 0.23 | 1/4 | 4mm | | 01-1205 |
| 33G | 0.2 | 1/4 | 4mm | | 01-1206 |

Remark: normal size in a above, we can also make special size as clients' request.

Product Conformance

Conform to ISO11608-2
In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PE,PP,SUS304 Stainless Steel Cannula,Silicone oil

