



## Biopsy Needles

### Brief description of the product

The Biopsy Needles are assembled by the inner needle (cutting needle) and an outer needle (cannula). The outer needle has graduation scale on it, and the inner needle has cutting edge and hopper. Biopsy needle is a medical device for kidney, liver, lung, breast, thyroid, prostate, pancreas, testes, uterus, ovaries, body surface and other organs. Biopsy needles tumor can be used for the sampling and cells drawing of cone tumors and unknown kind of tumors.

### Characteristic

- ◆ Omnidirectional 360° operation. Earth core positioning, puncture angle accurate to 1°. Target precision ±5mm. Cannula attached puncture needle is multifunction: biopsy, injection, body fluid collection, ablation single puncture, multitreatment. Assist minimally invasion.
- ◆ Biopsy: equipped with biopsy needles, single puncture, multi-sample collection, none resistance cutting, ensure good quality sample. medicine injection, prevent implantation metastasis and bleeding.

### Detailed specifications

SIZE		COLOR	LENGTH	REF
OD(mm)	GAUGE	COLOR & CODE	90 3 1/2"	
1.2	18G	pink		01-1401
1.4	17G	red purple		01-1402
1.6	16G			01-1403
1.8	15G	blue-grey		01-1404

Remark: normal size in above table, we can also make special size as clients' request.

### Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

### Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

### Main material

ABS,PP,PC, SUS304 Stainless Steel Cannula,Silicone oil

