




Scalp Vein Sets

Brief description of the product

The Scalp Vein Sets are assembled by Protective cap, needle tube, needle handle, tubing, needle hub. All materials used for the products meet medical requirement, sterilized by ETO,pyrogen-free. The Scalp Vein Sets intended to be used with Syringes or Infusion/Transfusion Sets for medicine infusion and blood collection.



Characteristic

- ◆ Normal type
- ◆ Soft wings
- ◆ Thin wall cannula, large ID, and smooth flow
- ◆ Safety Scalp Vein Sets, additional features
- ◆ Color coded wings by size for Cclear recognition
- ◆ Protective sheath on the tube
- ◆ Easy activation of safety sheath, prevent accidental needle sticks
- ◆ Needle hub: 1. PVC hub(slip with cap)
2. ABS hub(luer lock with cap)
- ◆ Packing: in bulk, PE bag, in blister.

Detailed specifications

SIZE OD (mm)	GAUGE	COLOR &CODE	LENGTH				ID		Bevel angle		REF
			1/2	5/8	3/4	1	RW W	T	Long bevel	Short bevel	
0.40	27G	Grey									01-201
0.45	26G	brown									01-202
0.50	25G	orange									01-203
0.55	24G	purple									01-204
0.60	23G	blue									01-205
0.70	22G	black									01-206
0.80	21G	green									01-207
0.90	20G	yellow									01-208
1.10	19G	cream									01-209
1.20	18G	pink									01-210

Remark: normal size in above table, we can also make special size as clients' request.

Product Conformance

Conform to ISO 7864
In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula,Silicone oil

