



Blood Collection(Multi-sample)Needles

Brief description of the product

The Blood-Collecting Needles are assembled by hub, needle tube, protect cap and rubber cap, All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free. The Blood-collecting needles intended to be used to blood or plasm collection.

Characteristic

1.Pen Type Blood Collection Needle: (Type I)

- ◆ Latex free
- ◆ Multi-sample needles permit several samples to be taken with a single puncture
- ◆ Sharp & smooth edges make penetration painless, easy connection to rubber stoppers

2.Blood Collection Needle with Scalp Vein Set (Type II) , Blood Collection Needle with Hypodermic Needle(Type III)

- ◆ Latex free
- ◆ Soft and transparent tube, can observe the vein blood flow clearly.
- ◆ Double wings make puncture more safe .
- ◆ Different color wings means different size.

Detailed specifications

SIZE		COLOR	LENGTH			ID		Bevel angle		REF
OD (mm)	GAUGE	COLOR & CODE	3/4	1	1 1/2	RW	TW	Long bevel	Short bevel	
0.60	23G	blue								01-401
0.70	22G	black								01-402
0.80	21G	green								01-403
0.90	20G	yellow								01-404
1.20	18G	pink								01-405

Remark: normal size in a above table, we can also make special size as clients' request.



SIZE		COLOR	LENGTH				ID		Bevel angle		REF
OD (mm)	GAUGE	COLOR & CODE	1/2	5/8	3/4	1	RW	TW	Long bevel	Short bevel	
0.55	24G	purple									01-406
0.60	23G	blue									01-407
0.70	22G	black									01-408
0.80	21G	green									01-409
0.90	20G	yellow									01-410
1.10	19G	cream									01-411
1.20	18G	pink									01-412

Remark: normal size in above table, we can also make special size as clients' request.



Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)



Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

ABS,rubber, SUS304 Stainless Steel Cannula,Silicone oil

