



Pre-Filled Flush Syringes

Brief description of the product

The Pre-Filled Flush Syringes are assembled by barrel, plunger, piston, syringe cap and Saline. All parts and material for this product meet medical requirements, after sterilized by moist heat. It's used to seal and wash the catheter end during the internal of different drug treatment.

Characteristic

- ◆ Three sizes comply different flush application
- ◆ BIIR piston, latex free
- ◆ Special design to prevent plunger slipping out.

Detailed specifications

TYPE	REF
3ml	02-1001
5ml	02-1002
10ml	02-1003

Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, Saline, BIIR, Silicone oil

