



NANCHANG KINDLY
MEDICAL TECHNOLOGY CO.,LTD.



Safety Needles


Brief description of the product

The Safety Needles are assembled by needle hub, protection arm, protection cap, needle tube. All materials used for the products meet medical requirement, sterilized by ETO, pyrogen-free. The Safety Needles are used together with syringes or infusion (blood) device, applied to human skin, subcutaneous, intramuscular, intravenous, etc. for drug (blood) injection or blood extraction, avoid cross-infection effectively.

Characteristic

- ◆ The design of safe protector cap can shield the cannula after use, prevent needle point from stabbing nurse, patient or other people, which can effectively reduce infection.

Instruction for use



Product Conformance

Conform to ISO 7864 and ISO 23908

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula, Silicone oil

