



NANCHANG KINDLY
MEDICAL TECHNOLOGY CO.,LTD.



Sterile Syringes

Brief description of the product

The Sterile Syringes (with/without needle) are assembled by barrel, plunger, piston, with/without Hypodermic needles, which can be classified to central nozzle type and side nozzle type according to the prick connector. It has three parts and two parts according to the structure. All parts and material for this product meet medical requirements, after sterilized by ETO, pyrogen free.

The Sterile Syringes used together with Hypodermic needles are intended to inject the drug for the patient or take out the blood for the patient.

Characteristic

- ◆ Different size, meet different requirements.
- ◆ Central luer slip/central luer lock/off-center luer slip.
- ◆ Special design to prevent plunger slipping out.
- ◆ Latex/latex free piston.
- ◆ W/needle, W/O needle.
- ◆ Custom needle size per specification.

Detailed specifications

SIZE	SLIP REF	LOCK REF	OFF-CENTER REF
1ml	02-101	02-112	-----
2ml	02-102	02-113	-----
3ml	02-103	02-114	-----
5ml	02-104	02-115	-----
10ml	02-105	02-116	02-123
20ml	02-106	02-117	02-124
25ml	02-107	02-118	02-125
30ml	02-108	02-119	02-126
35ml	02-109	02-120	02-127
50ml	02-110	02-121	02-128
100ml	02-111	02-122	-----

Product Conformance

Conform to ISO 7886-1 and ISO 7864

In compliance with European Medical Device Directive 93/42/EEC(CE Class: W/needle: IIa;W/O needle: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PP, SUS304 Stainless Steel Cannula,Isoprene rubber,Silicone oil

