



NANCHANG KINDLY
MEDICAL TECHNOLOGY CO.,LTD.



Feeding Tubes

Brief description of the product

The Feeding Tubes are assembled by tube, basic pedestal and pedestal stopple. All parts and material for this product meet medical requirements, after sterilized by ETO. It is mainly used for patient who can't take any food to get the necessary nutrition.

Characteristic

- ◆ Medical grade PVC
- ◆ Appropriate degree of hardness provides kink resistance
- ◆ Soft rounded closed tip
- ◆ Two lateral eyes with smooth edges
- ◆ Color coded connector

Detailed specifications

SIZE	REF
4Fr	05-701
6Fr	05-702
8Fr	05-703
10Fr	05-704

SIZE	REF
4Fr	05-705
6Fr	05-706
8Fr	05-707
10Fr	05-708

Product Conformance

Conform to EN 1615

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PVC



① 701-704

② 705-708