



Suction Catheter Kits

Brief description of the product

The contents in the Suction Catheter kits includes 1 each catheter, 1 each vinyl glove, 1 each pop-up cup.

The Suction Catheter is made up of tube and connector with thumb control. It is sterile. It's mainly used for removal of secretion from trachea and bronchial tube to avoiding choke and death.

Characteristic

- ◆ Medical grade PVC tubing
- ◆ ABS valve
- ◆ Pop-up solution cup
- ◆ 1 or 2 vinyl glove included
- ◆ Sterile

Detailed specifications

SIZE	LENGH	GRADUATION	REF
6Fr	435mm	with	07-201
8Fr	435mm	with	07-202
10Fr	585mm	with	07-203
12Fr	585mm	no	07-204
14Fr	585mm	no	07-205
16Fr	585mm	no	07-206
18Fr	585mm	no	07-207

Product Conformance

Conform to EN 1733

In compliance with European Medical Device Directive 93/42/EEC(CE Class: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

PVC,ABS,Vinyl Glove,Paper

