



NANCHANG KINDLY
MEDICAL TECHNOLOGY CO.,LTD.



Withdrawal and Injection Spikes

Brief description of the product

The Withdrawal and Injection Spikes are assembled by protective cap, spike, filter, adapter. used to draw up medicine liquid from container, with features such as filtering, prevent counter current.

Characteristic

- ◆ Medical grade material
- ◆ Precision filtration
- ◆ Sterile

Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

ABS,PES,PP

