

Aesthetic Products

Sterile Aesthetic Cannula

Blunt end, side hole, with scale

Sterile Syringe

Luer lock, slender diameter, high injection pressure, suitable for injection of high viscosity filler.

The finger is designed to fit the large hand peg, and the annular push rod is easy to operate with one hand, which can effectively control the injection speed. High permeability medical materials with good drug compatibility.

Hypodermic Needle

Microneedle, Size: 30G-34G

Certificate and Quality Assurance:

CE, ISO13485, FDA

Main material

PP / PC, Isoprene Rubber, Silicone oil, SUS304 Stainless Steel Cannula

Intend use

These products are intended for hypodermic injection of body fluids or tissues, drugs or medical devices approved for human use.



Three Ring Sterile syringes

Main material

PP, Isoprene Rubber, Silicone oil, ABS

Certificate and Quality Assurance:

CE, ISO13485

Intend use

Sterile syringes are intended to inject drug for patients.

And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

Luer lock, suitable for injection of high viscosity filler.

The finger is designed to fit the large hand peg, and the annula push rod is easy to operate with one hand, which can effectively control the injection speed.

