

Disposable Bone Marrow Infusion Sets for Single Use

Needle set, needle tip protector, needle set protector, sealed dialyzed paper .

Certification and Quality Assurance

Conform to ISO11608-2

In compliance with European Medical Device Directive 93/42/EEC(CE Class: IIa)

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main Material

PE, PP, SUS304 Stainless Steel Cannula, Silicone Oil

Description

- 29-33G, needle length 4mm-12mm, thin wall/regular wall
- Double point system
- Sterile, non-toxic. non-pyrogenic
- Safety design and easy for using
- High compatibility and cover almost all branches

Needle Length	4mm-12mm
Needle Size	29-33G

