

Disposable Multi-Needles (9 Pin)

Product features:

- Manufactured from high-quality austenitic stainless steel for strength and durability
- Thin-wall needle tube design offers a larger inner diameter, allowing higher flow rates
- Precision-engineered needle tip ensures exceptional sharpness for fast insertion, reduced pain, and minimal tissue trauma
- 6:100 conical connector complies with international standards, ensuring precise sizing and excellent compatibility with other medical devices



Structure and composton

Multi-needle is composed of needle hub, needle tube and protective cap.

Intended use

Multi-needle is intended to inject fluids into the human body.

Main Material

PP, SUS304 Stainless Steel Cannula, Silicone Oil

Certification and Quality Assurance

Manufacturing process is in compliance with ISO 13485 Quality System.

OD	Gauge	General Specifications
0.2	33G	33G
0.23	32G	32G
0.25	31G	31G
0.30	30G	30G

Note: the specification and length can be customized according to the customers’ requirements