



Enteral Product Catalogue



Shanghai Kindly(KDL) Group
Nanchang Kindly Medical Technology Co.,Ltd.

Add: F/7,Tower A,No.77 North Aixihu Road, Nanchang High-Tec Industry Development
Zone, Nanchang, Jiangxi, China

Post Code: 330096

Tel: +86-0791-86861216

Mob: +86-13922143599

Web: www.kdInc.com

E-mail: Melvin.lu@kdInc.com


Rev . : 20191130

To address the risk of wrong route administration ISO 80369-3 medical device connectors are being progressively introduced worldwide, replacing conventional Luer connectors. Sampling and delivery devices intended for enteral applications and for major regional enteral will in future utilise ISO 80369-3 connectors.

Enteral connectors devices will for easy identification be coloured purple.

Product Catalogue

Enteral Syringe	01/02
Enteral Feeding Spike Set	03/04
Enteral Feeding & Drainage Tube	05/06
Enteral Feeding Bag	07/08
Oral Syringe	09/10
Bottle Adaptor	11/12



Enteral Syringe

Structure and composition:

The Enteral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.

Material:

PP, Isoprene rubber, Silicone oil

The Enteral syringe is used for delivering drug or food to oral or enteral.

Specifications:

1ml, 3ml, 5ml, 10ml, 20ml, 60ml.

Coloured code purple

Enteral syringes with female EN
Fit connection (ISO 80369-3)



Enteral Single Use Syringes (Isoprene rubber piston)



Enteral Home Use Syringes (Silicone gasket)

DEHP
FREE

LATEX
FREE

Enteral Feeding Spike Set

Enteral feeding spike set, male EN Fit connector (ISO 80369-3) And female EN Fit connector.
It is suitable for enteral feeding patients. Both inpatients and non inpatients can use this product for feeding,
For enteral feeding only.


Product Features:

**DEHP
FREE**

CLAMP

**LATEX
FREE**

**ENFit
COMPLIANT**



Good compatibility, the interface can be quickly linked with the bottled preparations in bags and large and small bottles.




Enteral Feeding & Drainage Tube


The EN-fit Feeding & Drainage Tube made from PVC, has clear numbered cm markings and is fitted with the new ENFit (ISO 80369-3) connector.

It is suitable for enteral feeding patients. Both inpatients and non inpatients can use this product for feeding.
For enteral feeding only

Clear numbered cm markings




Extension set/syringe driver connection set, male EN Fit connector (ISO 80369-3) And female EN Fit connector.



Enteral Feeding Bag

Enteral Feeding Bag is a durable enteral feeding bag that comes with attached Enteral feeding spike set consisting of flexible drip chamber pump set or gravity set, built-in hangers and a large top fill opening with leak-proof cap. The bags are made from soft EVA material.



Extension set/syringe driver connection set, male EN Fit connector (ISO 80369-3) And female EN Fit connector.



Oral Syringe

Structure and composition:

The Oral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.

Material:

PP, Silicone gasket, Silicone oil

Oral syringe for single use are intended for oral administration of drugs.

Specifications:

1ml, 3ml, 5ml, 10ml, 20ml.



Bottle Adaptor

Material:
PE, TPE

Suitable with Enteral syringe



Suitable with Oral syringe

