


Huber Needle

Structure and composition:

The huber needles are assembled by protective cap, needle hub, needle tube, tubing, injection site, robert clamp and other parts. All materials meet medical requirement, sterilized by ETO, pyrogen-free.

Material:

PC, PP, PVC, SUS304 Stainless Steel Cannula, Silicone oil



Remark: The type can make as clients' request.

Intended use:

The huber needle is intended to be used to feed medicine to the device implanted in patient hypoderm.

Detailed specifications

SIZE		COLOR	LENGTH								ID		Bevel	angle
OD (mm)	GAUGE	COLOR & CODE	1/3"	1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"	RW	TW	Long bevel	Short bevel
0.70	22G	black												
0.80	21G	green												
0.90	20G	yellow												
1.10	19G	cream												

Remark: normal size in above table, we can also make special size as clients' request.

