

Disposable Infusion Sets (Plastic Spike)

The product consists of a corkscrew puncture device with protective sleeve, plastic puncture device, air filter or intake device, dropper, tubing, stopcock, tee connector, injection components, flow regulator, medication filter, external cone connector with protective sleeve, and intravenous infusion needle.

Product Use

For intravenous gravity infusion of medication to patients.

Product features:

- Complete specifications to meet the patient's infusion requirements;
- The cone connector is designed according to the 6:100 standard, which is suitable for medical equipment.

Product Specifications:

- Inlet type is suitable for infusion bottles;
- Non-inlet type is suitable for infusion bags.

Certification and Quality Assurance

Manufactured in full compliance with ISO13485 quality systems.

Needle specifications: 0.4, 0.45, 0.5, 0.55, 0.6, 0.7, 0.8, 0.9, 1.1, 1.2mm.

