



Cosmetic needle

Size	Lenth(mm)				Color mark
18G				70	● pink
22G		40	50	70	● black
23G			50		● blue
25G		40	50		● orange
27G		40			● grey
30G	25				● yellow
31G	25				○ white

Micro needle

Size	Lenth(mm)				Color mark
30G	1.5	2.5	4	> 4mm	● yellow
31G	1.5	2.5	4	> 4mm	○ white
32G	1.5	2.5	4	> 4mm	● green
33G	1.5	2.5	4	> 4mm	● black
34G	1.5	2.5	4	> 4mm	● orange

Remark: normal size in above table, we can also make special size as clients' request.

Shanghai Kindly Group

Nanchang Kindly Medical Technology Co.,Ltd.

Add: F/7,Tower A,No.77 North Aixihu Road, Nanchang High-Tec Industry Development Zone, Nanchang, Jiangxi. China

Post Code: 330096


Tel: +86-0791-86861216

Mob: +86-13922143599

Web: www.kdInc.com

E-mail: Melvin.lu@kdInc.com

V1.0-2024




TO WIN UNIVERSAL CONFIDENCE WITH KDL QUALITY AND REPUTATION




COSMETIC NEEDLE



PRODUCT STRUCTURE



PRODUCT ADVANTAGE



- 1. Effectively avoid the risk of tissue trauma caused by direct filling of traditional sharp needles, and prevent the occurrence of other complications caused by sodium hyaluronate entering blood vessels and causing embolism.
- 2. Effective pain relief; The rounded and blunt design of the needle tip avoids multiple punctures when sliding between tissues and injuring blood vessels and nerves.
- 3. Full facial injection possible. Especially in sensitive areas (around the eyes, tip of the nose, around the temples), blunt needles have a unique advantage.

PRODUCT ADVANTAGE



MULTIPLE PACKAGING

