

Sterile Syringes for Single Use (3-Parts)

Structure and Composition:

The Syringes are assembled by Barrel, Plunger, Piston with/without Hypodermic Needles.
All parts and material for this product meet medical requirements. Sterilized by EO.

Main Material:

PP, Isoprene Rubber, Silicone Oil

Intend Use:

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

Certificate and Quality Assurance:

CE, FDA, ISO13485



Sterile Syringes for Single Use (Luer Slip)

Specification:

Luer Slip, Luer Lock

1ml, 2ml, 2.5ml, 3ml, 5ml, 10ml, 20ml, 25ml, 30ml, 35ml, 50ml, 60ml, 100ml

Refer Code:

| Specification | Luer Lock | Luer Slip | Eccentric Tip |
|---------------|-----------|-----------|---------------|
| 1ml | SR11101 | SR11201 | |
| 2ml | SR11102 | SR11202 | |
| 2.5ml | SR1112.5 | SR1122.5 | |
| 3ml | SR11103 | SR11203 | |
| 5ml | SR11105 | SR11205 | |
| 10ml | SR11110 | SR11210 | SR11310 |
| 20ml | SR11120 | SR11220 | SR11320 |
| 30ml | SR11130 | SR11230 | SR11330 |
| 35ml | SR11135 | SR11235 | SR11335 |
| 50ml | SR11150 | SR11250 | SR11350 |
| 60ml | SR11160 | SR11260 | SR11360 |
| 100ml | SR111100 | SR112100 | SR113100 |



Suitable to be connected with power-driven syringe pumps



Sterile Syringes for Single Use (Luer Lock)