

## Low Dead Space Needle

### Description:

- This needle incorporates a Low Dead Space (LDS) hub which is design to reduce the effective space left behind in the needle hub close to zero.
- This new patented LDS hub design can result in savings of up to 0.08ml which leads to significant cost savings when injecting high cost drugs.
- Injecting high cost drugs
- Injecting Anti-VEGF for the treatment of macular edema
- Ophthalmic treatments
- Vaccination
- Cosmetic injections, especially of the face, neck, hands
- Local Anesthesia

### Material

- PVC/stainless steel needle.

### Certificate

- ISO,CE

