



Oral Rinsing Needles

The product, a disposable, non-sterile oral irrigation system, consists of a syringe, needle holder, and an optional positioning device. It requires sterilization prior to use as per the instructions for use.

Medical institutions use it for removing debris or foreign objects in the mouth during oral treatment.

Certification and Quality Assurance

In compliance with Medical Devices Directive 93/42/EEC(Class IIa)

Manufacturing process is in compliance with ISO 13485

Main Material

PP, SUS304

Needle Size		Specification	
Gauge	mm	Tip Type	Wall Type
31G	0.25	Round,Flat,or Beveled	Regular Wall (RW), Thin Wall (TW)
30G	0.3		
29G	0.33		
28G	0.36		
27G	0.4		
26G	0.45		
25G	0.5		