



# Bone Marrow Biopsy Needles

**SHANGHAI KINDLY(KDL) GROUP**  
NANCHANG KINDLY MEDICAL TECHNOLOGY CO.,LTD

F: +86-791 86861216 W: [www.kdInc.com](http://www.kdInc.com)  
T: +86-791 86861216

The handle includes a sharpened luer-lock which gives a safer Attachment to the syringe. Aspiration can be performed at any Step of the procedure


Digital scale can provide accurate puncture depth

### Multiple edge

The edge is strong and sharp, that allows a quick and less invasive penetration into the cortical

Special trocar needle and sharp tubing ensure minimal invasion, 5 side holes at front end of tubing, multi-way collect bone marrow, more efficient and less pain.

Pick up needle is provided with a rotation coupling security system which limits the risk of movement of the stylet during the insertion.



## Structure and composition:

The Bone marrow biopsy needle consist of needle hub, needle tube, lining core, handle and Protective cap.

## Material:

ABS,PP PC,SUS304 Stainless Steel Cannula,Silicone oil

## Intended use:

The Bone marrow biopsy needle is used for bone marrow aspiration.



## Specification

| SIZE   |       | LENGTH     |            |             |           |
|--------|-------|------------|------------|-------------|-----------|
| OD(mm) | GAUGE | 70<br>2 ¾" | 90<br>3 ½" | 110<br>3 ¼" | 130<br>5" |
| 1.8    | 15G   | ✓          |            | ✓           |           |
| 2.5    | 13G   | ✓          |            | ✓           | ✓         |
| 3.0    | 11G   |            | ✓          |             |           |

Remark: normal size in above table, we can also make special size as clients' request.