



Sterile Syringes for Insulin for Single Use (Retractable)

Structure and Composition:

The Syringes are assembled by Sliding Sleeve, Needle Protective Cap, Needle Tube, Barrel, Plunge and Piston. All parts and material for this product meet medical requirements. Sterilized by EO.

Main Material:

PP, Isoprene Rubber, Silicone Oil, SUS304 Stainless Steel Cannula

Specification:

1ml, 0.5ml, 0.3ml
U-40, U-100

Refer Code:

Specification	U-40	U-100
0.3ml	INS211	/
0.5ml	INS221	INS222
1ml	INS231	INS232

Needle Size:

27G-31G



Certificate and Quality Assurance:

CE, FDA, ISO13485

Intend Use:

The sterile syringes for insulin for single use intended to be used to administer an injection of insulin to a patient subcutaneously.