

# Sterile Syringes for Single Use (2-Parts)

## Structure and Composition:

The Syringes are assembled by Barrel, Plunge, with/without Hypodermic Needles.  
All parts and material for this product meet medical requirements. Sterilized by EO.

## Main Material:

PP, Silicone Oil, SUS304 Stainless Steel Cannula

## Specification:

Luer Slip

1ml, 2ml, 3ml, 5ml, 10ml, 20ml

## Refer Code:

Specification	Luer Lock
1ml	SR21201
2ml	SR21202
3ml	SR21203
5ml	SR21205
10ml	SR21210
20ml	SR21220



## Certificate and Quality Assurance:

CE, ISO13485

## Intend Use:

Sterile syringes with needle are intended to inject drug for patient. And syringes are intended to use immediately after filling and are not intended to contain the medicament for extended periods of time.