

# Sterile Syringes for Single Use (Colorful Plungers)

## Structure and Composition:

The Syringes are assembled by Barrel, Plunge, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

## Main Material:

PP, Isoprene Rubber, Silicone Oil

## Specification:

Luer Slip, Luer Lock

1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml

## Certificate and Quality Assurance:

CE, ISO13485

## Refer Code:

Specification	Luer Lock	Luer Slip	Eccentric Tip
1ml	SR11101	SR11201	
2ml	SR11102	SR11202	
2.5ml	SR1112.5	SR1122.5	
3ml	SR11103	SR11203	
5ml	SR11105	SR11205	
10ml	SR11110	SR11210	SR11310
20ml	SR11120	SR11220	SR11320
30ml	SR11130	SR11230	SR11330
35ml	SR11135	SR11235	SR11335
50ml	SR11150	SR11250	SR11350
60ml	SR11160	SR11260	SR11360
100ml	SR111100	SR112100	SR113100



## Intend Use:

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.