



Sterile Syringes for Single Use (Fixed Needles)

Structure and Composition:

The Syringes are assembled by Barrel, Plunger, Piston, Fixed Needles.
All parts and material for this product meet medical requirements. Sterilized by EO.

Certificate and Quality Assurance:

CE, FDA, ISO13485

Main Material:

PP, Isoprene Rubber, Silicone Oil,
SUS304 Stainless Steel Cannula

Needle Size:

27G-31G

Specification:

1ml

Refer Code:

1ml SR11401



Fixed Needle

Intend Use:

Sterile syringes are intended to inject drug for patients.
And syringes are intended for use immediately after filling
and are not intended to contain the medicament for
extended periods of time.