

Sterile Syringes for Single Use (PC)

Structure and Composition:

The Syringes are assembled by Barrel, Plunge, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

Main Material:

PC, Isoprene Rubber, Silicone Oil

Specification:

Luer Lock

1ml, 3ml, 5ml, 10ml, 20ml, 30ml

Refer Code:

Specification	Luer Lock
1ml	SR12101
3ml	SR12103
5ml	SR12105
10ml	SR12110
20ml	SR12120
30ml	SR12130



Certificate and Quality Assurance:

CE, ISO13485, FDA

Intend Use:

Sterile syringes with needle are intended to inject drug for patient. And syringes are intended to use immediately after filling and are not intended to contain the medicament for extended periods of time.