

Sterile Syringes for Veterinary Use

Structure and Composition:

The Syringes are assembled by Barrel, Plunger, Piston, Protective Cap, with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

Main Material:

PP, Isoprene Rubber, Silicone Oil, sus304 Stainless Steel Cannula

Refer Code:

Specification	Luer Lock
3ml	VSR103
5ml	VSR105
10ml	VSR110
20ml	VSR120
30ml	VSR130
60ml	VSR160

Specification:

With Needles: 3ml, 5ml, 10ml

Without Needles: 3ml, 5ml, 10ml,
20ml, 30ml, 60ml



Veterinary Needle

Specification :

Plastic hub

30G	Effective Length $13 \leq L \leq 50\text{mm}$
27G	
26G	
25G	
24G	
23G	
22G	
21G	
20G	

Aluminum Hub

18G	Effective Length $13 \leq L \leq 50\text{mm}$
16G	
14G	

