



Transfer Spikes

Brief description of the product

The Transfer Spikes are assembled by plastic spikes, protective caps. All parts and material for this product meet medical requirements, after sterilized by ETO, pyrogen free.

Transfer Spikes are intended for transfer or mixing of sterile fluids/drugs in closed systems, for solubilization of lyophilisation or suspension of erythrocytes.

Characteristic

- ◆ Umbilical cord clamp
- ◆ Individually packaged
- ◆ Sterile

Detailed specifications

TYPE	REF
Green	08-101
Blue	08-102

Product Conformance

In compliance with European Medical Device Directive 93/42/EEC(CE Class: I)

Quality Assurance

Manufacturing process is in compliance with ISO 13485 and ISO9001 Quality System.

Main material

ABS, PP

