

Vacuum Urine Container

Description:

- The optimized threading of the urine cup makes it easy to open and close.
- Has a leak-resistant screw cap to minimize risk of contamination and exposure to specimen.
- It allows urine to be transferred into one or more vacuum tube without exposure to the sample.
- Reducing leakage,spills and contamination.

Material

- PP/stainless steel needle.

Certificate

- ISO,CE



Item No	Description of Goods	Quantity	Packing	Measurement(cm)
GDKK1027	Urine transfer device for urine tube,100ml-120ml	250	1pc*250bags	60*31*44
GDK0100UA	Vacuum Urine Tube,Yellow cap,no Additive,16*100mm,10ml	1200	100pc*12racks	51*40*24
GDK-H4	Urine Transfer Straw,White,PE	5000	50pcs/bag*100bags	38*35*17

