

## Description

- The urine cup features optimized threading for easy opening and closing, and is equipped with a leak-resistant screw cap to minimize the risk of contamination and specimen exposure. The integrated transfer device ensures precise measurements, enhancing accuracy during testing.



## Vacuum Urine Cup

This product is designed for the collection of liquid samples during the inspection process and is compatible with vacuum tubes. It is managed in accordance with non-medical device regulations.

- Made from durable hard plastic, the cup maintains its shape under normal usage conditions. It is marked with clear scales for easy observation of liquid collection levels, facilitating better storage and handling.
- The special sealing design, when used in conjunction with vacuum tubes, ensures that liquids are securely contained with no exposure to the outside environment, preventing odors from escaping and improving the laboratory working conditions.