


we desire every patients
with diabetes
can be healthy



Insulin pen needle contrast before and after using



Before using After once use After twice use After sixth use

Type Specification

Type	Specification
33G NEW	0.20x4mm
32G NEW	0.23x4mm
31G	0.25x5mm
31G	0.25x8mm
30G	0.30x5mm
30G	0.30x8mm
29G	0.33x12mm

Shanghai Kindly(KDL) Group
Nanchang Kindly Medical Technology Co.,Ltd
Add:F/7,Tower A,No.77 North Aixihu Road,Nanchang High-Tec Industry
Development Zone,Nanchang,Jiangxi.China.
Tel: +86-0791-86861216
Web: www.kdInc.com
E-mail: inquiry@kdInc.com



Disposable
Insulin Pen Needle

Insulin Pen Needle Features

- Commonality** Suit any kind of insulin pen and easy for medical workers and patients to use.
- All specifications** The product is complete in specification and can satisfy different injection requirement of patients.
- Comfortableness** Reduce painful feeling of puncture and make patients feel comfortable during injection.
- Professional** KDL is one of the largest medical needle manufacturers in China.

Product Features



Safety design and easy for using

The strength design in needle set make it more convenient and safe to use


High compatibility and cover almost all branches

The product has high compatibility. The standard design of internal threaded ensure it can stably match all branches insulin pen.

Best penetration make the injection more comfortable

KDL insulin pen needle has a series of advanced and mature needle tip technique and silication technique. The sharp needle tip make puncture more easily and let patients feel more comfortable during injection.

High Matching Degree



The screw thread design of KDL insulin pen needle meets the standards regulated by ISO11608-2, which is compatible with other standard ones domestically produced or imported.

More Comfortable, More Safety

- Sharpening Technique**
International sharpening technique make KDL needle tip has excellent penetration.
- Silication Technique**
The advanced silication technique can lower resistance of needle tip and cannula. Reduce painful feeling of injection and make patients feel more comfortable.
- Thin Wall Stainless Material**
KDL insulin pen needle use thin wall medical stainless material whose inner diameter is larger than the stainless needle with same outer diameter. It ensure insulin injection more fluently and organism can absorb the insulin rapidly.



Needle Structure Parsing

Sterile injection guarantee with double protector



33G Superfine Short Insulin Pen Needle

4mm needle is 50% shorter than traditional 8mm needle


NEW

Reduce the risk of injecting to muscle layer directly and make the injection more safely.

Don't need to pinch the skin part during injection which make it more convenient and easy in operation.

33G insulin pen needle will create more comfortable injection feeling to patient.

33G, a shorter and finer insulin pen needle!!



Don't need to pinch the skin when use 4mm needle

Contrast of 8mm and 4mm needle used in few subcutaneous tissue area